

ASTORIA - PACIFIC



Astoria-Pacific[™]
INTERNATIONAL

Automated analysis made even easier



Easier And Less Expensive To Use

- Less time to train personnel
- Easier to maintain
- Less reagent usage
- Automatically makes calibrants from a stock standard
- Automatically dilutes out-of-range samples
- Automatically detects low levels of reagents or samples



rAPID Software

This data acquisition software for discrete analysis has been designed with the operator in mind. The built-in tracking features keep the operator visibly informed of the status of each reagent, sample and reaction well.

Sample tables can be imported from LIMS; reports can be easily customized and exported to LIMS as ASCII files or spreadsheet applications.

Also, rAPID works within the WINDOWS® environment; therefore its useful features, like auto spikes, duplicates, and QC charts, are only a click away.

Analytical Applications

The Astoria Discrete is a flexible platform, offering tests that satisfy EPA, AOAC and other regulatory standards. Enzymatic methods, such as those in the wine and beverage industry, can also be performed on the Astoria Discrete.

The development of new methods through rAPID is a simple process, thereby increasing the future versatility of the Astoria Discrete in your laboratory

Reagent Kits

Kits for specific methods are available, manufactured to cGMP quality standards.

Features

- Dual syringes offer high precision sampling and dispensing
- Rack for 27 reagents: 15 x 15-mL bottles, 12 x 30-mL bottles
- 60 or 96 position sample racks available
- 96-place reaction plate, with temperature control
- Automatic QC of reaction wells
- Four built in detectors
- Automatic dilution and calibration preparation built in
- Separate wash stations for probe and reaction wells
- Compact, lightweight design
- Low reagent and sample usage

For more information on the Astoria Discrete, please contact us directly or through your local representative.