

NEW for Wine and Food Laboratories: Astoria•Pacific manual test kits

Test your samples immediately!

Reduce outside testing!

Test more frequently!

Multiple tests to choose from!

(Acetic Acid, Glucose + Fructose, L-Malic Acid and more!)



Introducing the Awareness Stat Fax 4500:

- Accommodates 12mm tubes or cuvettes.
- Interactive touch-screen LCD
- Read cell is temperature controlled to 37°C
- Bichromatic reading capability
- Standard model includes six filters - 340, 405, 505, 545, 580, and 630nm (or custom)
- On board curve fitting calculations with step-by-step user prompting.
- Single point calibration by standard or factor, multipoint calibration, kinetic rate by standard or factor (batch or singly).
- Large non-volatile memory stores 120 user entered tests.

